

PRKAG3 Antibody



PACO43612

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

PRKAG3

Uniprot

Q9UGI9

Synonyms:

5'-AMP-activated protein kinase subunit gamma-3 (AMPK gamma3) (AMPK subunit gamma-3), PRKAG3, AMPKG3


Immunogen:

Recombinant Human 5'-AMP-activated protein kinase subunit gamma-3 protein (1-210AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: 5-AMP-activated protein kinase subunit gamma-3 antibody at 5 μ g/ml. Lane 1: Hela whole cell lysate. Lane 2: HepG2 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 55, 52 kDa. Observed band size: 55 kDa.