

ANXA7 Antibody



PACO43673

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Gene ID:

ANXA7

Uniprot

P20073

Synonyms:

Annexin A7 (Annexin VII) (Annexin-7) (Synexin), ANXA7, ANX7 SNX


Immunogen:

Recombinant Human Annexin A7 protein (259-488AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: ANXA7 antibody at 4 μ g/ml. Lane 1: Jurkat whole cell lysate. Lane 2: LO2 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 53, 51 kDa. Observed band size: 53 kDa.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO43673 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43673 at dilution of 1:100.