

APOL1 Antibody



PACO43677

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May play a role in lipid exchange and transport throughout the body. May participate in reverse cholesterol transport from peripheral cells to the liver.

Gene ID:

APOL1

Uniprot

O14791

Synonyms:

Apolipoprotein L1 (Apolipoprotein L) (Apo-L) (ApoL) (Apolipoprotein L-I) (ApoL-I), APOL1, APOL


Immunogen:

Recombinant Human Apolipoprotein L1 protein (179-398AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43677 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43677 at dilution of 1:100.