

BMX Antibody



PACO43686

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

BMX

Uniprot

P51813

Synonyms:

Cytoplasmic tyrosine-protein kinase BMX (EC 2.7.10.2) (Bone marrow tyrosine kinase gene in chromosome X protein) (Epithelial and endothelial tyrosine kinase) (ETK) (NTK38), BMX


Immunogen:

Recombinant Human Cytoplasmic tyrosine-protein kinase BMX protein (65-280AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: Cytoplasmic tyrosine-protein kinase BMX antibody at 4 μ g/ml. Lane 1: A549 whole cell lysate. Lane 2: Mouse heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 78 kDa. Observed band size: 78 kDa.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43686 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon tissue using PACO43686 at dilution of 1:100.