

CDH16 Antibody



PACO43699

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Gene ID:

CDH16

Uniprot

O75309

Synonyms:

Cadherin-16 (Kidney-specific cadherin) (Ksp-cadherin), CDH16


Immunogen:

Recombinant Human Cadherin-16 protein (19-230AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43699 at dilution of 1:100.