
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in lifespan determination in ubiquinone-independent manner. Involved in ubiquinone biosynthesis. Potential central metabolic regulator.

Gene ID:

COQ7

Uniprot

Q99807

Synonyms:

5-demethoxyubiquinone hydroxylase, mitochondrial (DMQ hydroxylase) (EC 1.14.13. -) (Timing protein clk-1 homolog) (Ubiquinone biosynthesis monooxygenase COQ7), COQ7


Immunogen:

Recombinant Human 5-demethoxyubiquinone hydroxylase, mitochondrial protein (37-217AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43712 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO43712 at dilution of 1:100.