

AGGF1 Antibody



PACO43718

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Promotes angiogenesis and the proliferation of endothelial cells. Able to bind to endothelial cells and promote cell proliferation, suggesting that it may act in an autocrine fashion.

Gene ID:

AGGF1

Uniprot

Q8N302

Synonyms:

Angiogenic factor with G patch and FHA domains 1 (Angiogenic factor VG5Q) (hVG5Q) (G patch domain-containing protein 7) (Vasculogenesis gene on 5q protein), AGGF1, GPATC7 GPATCH7 VG5Q


Immunogen:

Recombinant Human Angiogenic factor with G patch and FHA domains 1 protein (1-70AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43718 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO43718 at dilution of 1:100.