

FKBP2 Antibody



PACO43741

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:200-1:1000,
IHC:1:20-1:200

Protein Background:

PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

Gene ID:

FKBP2

Uniprot

P26885

Synonyms:

Peptidyl-prolyl cis-trans isomerase FKBP2 (PPlase FKBP2) (EC 5.2.1.8) (13 kDa FK506-binding protein) (13 kDa FKBP) (FKBP-13) (FK506-binding protein 2) (FKBP-2) (Immunophilin FKBP13) (Rotamase), FKBP2, FKBP13


Immunogen:

Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP2 protein (22-142AA).

Storage:


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot.

All lanes: Peptidyl-prolyl cis-trans isomerase FKBP2 antibody at 10µg/ml. Lane 1: Mouse thymus tissue. Lane 2: Rat brain tissue. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 16 kDa. Observed band size: 16 kDa.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43741 at dilution of 1:100.