

# LRG1 Antibody



PACO43751

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:****Gene ID:**

LRG1

**Uniprot**

P02750

**Synonyms:**

Leucine-rich alpha-2-glycoprotein (LRG), LRG1, LRG

**Immunogen:**


Recombinant Human Leucine-rich alpha-2-glycoprotein protein (36-347AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human lung cancer using PACO43751 at dilution of 1:100.