

MEP1A Antibody



PACO43752

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

MEP1A

Uniprot

Q16819

Synonyms:

Mepirin A subunit alpha (EC 3.4.24.18) (Endopeptidase-2) (N-benzoyl-L-tyrosyl-P-amino-benzoic acid, hydrolase subunit alpha) (PABA peptide hydrolase) (PPH alpha), MEP1A


Immunogen:

Recombinant Human Mepirin A subunit α protein (440-700AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43752 at dilution of 1:100.