

TGM2 Antibody



PACO43769

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

Gene ID:

TGM2

Uniprot

P21980

Synonyms:

Protein-glutamine gamma-glutamyltransferase 2 (EC 2.3.2.13) (Tissue transglutaminase) (Transglutaminase C) (TG(C)) (TGC) (TGase C) (Transglutaminase H) (TGase H) (Transglutaminase-2) (TGase-2), TGM2


Immunogen:

Recombinant Human Protein-glutamine gamma-glutamyltransferase 2 protein (438-687AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43769 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43769 at dilution of 1:100.