

AMBRA1 Antibody



PACO43782

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulates autophagy and development of the nervous system. Involved in autophagy in controlling protein turnover during neuronal development, and in regulating normal cell survival and proliferation.

Gene ID:

AMBRA1

Uniprot

Q9C0C7

Synonyms:

Activating molecule in BECN1-regulated autophagy protein 1, AMBRA1, KIAA1736

Immunogen:

Recombinant Human Activating molecule in BECN1-regulated autophagy protein 1 protein (1-240AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A