
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Induces activation of NF-kappa-B and stimulates production of the proinflammatory chemokine IL-8. Proinflammatory cytokine favoring Th2-type immune responses.

Gene ID:

IL25

Uniprot

Q9H293

Synonyms:

Interleukin-25 (IL-25) (Interleukin-17E) (IL-17E), IL25, IL17E


Immunogen:

Recombinant Human Interleukin-25 protein (33-177AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO43858 at dilution of 1:100.