

TOMM70 Antibody



PACO43865

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Receptor that accelerates the import of all mitochondrial precursor proteins.

Gene ID:

TOMM70

Uniprot

O94826

Synonyms:

Mitochondrial import receptor subunit TOM70 (Mitochondrial precursor proteins import receptor) (Translocase of outer membrane 70 kDa subunit) (Translocase of outer mitochondrial membrane protein 70), TOMM70, KIAA0719 TOM70 TOMM70A


Immunogen:

Recombinant Human Mitochondrial import receptor subunit TOM70 protein (60-340AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: TOMM70A antibody at 6 μ g/ml + MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 67 kDa. Observed band size: 67 kDa.