

MEF2C Antibody



PACO43903

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

MEF2C

Uniprot

Q06413

Synonyms:

Myocyte-specific enhancer factor 2C (Myocyte enhancer factor 2C), MEF2C


Immunogen:

Recombinant Human Myocyte-specific enhancer factor 2C protein (170-380AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO43903 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43903 at dilution of 1:100.