

PACO43909

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

Gene ID:

IFNA4

Uniprot

P05014

Synonyms:

Interferon alpha-4 (IFN-alpha-4) (Interferon alpha-4B) (Interferon alpha-76) (Interferon alpha-M1), IFNA4


Immunogen:

Recombinant Human Interferon alpha-4 protein (24-189AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human weasand tissue using PACO43909 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43909 at dilution of 1:100.