

SSX5 Antibody



PACO43924

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Could act as a modulator of transcription.

Gene ID:

SSX5

Uniprot

O60225

Synonyms:

Protein SSX5, SSX5


Immunogen:

Recombinant Human Protein SSX5 protein (1-229AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43924 at dilution of 1:100.