

CST6 Antibody



PACO43948

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Shows moderate inhibition of cathepsin B but is not active against cathepsin C.

Gene ID:

CST6

Uniprot

Q15828

Synonyms:

Cystatin-M (Cystatin-6) (Cystatin-E), CST6


Immunogen:

Recombinant Human Cystatin-M protein (29-149AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human skin tissue using PACO43948 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO43948 at dilution of 1:100.