

FOSL2 Antibody



PACO43965

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Controls osteoclast survival and size. As a dimer with JUN, activates LIF transcription. Activates CEBPB transcription in PGE2-activated osteoblasts.

Gene ID:

FOSL2

Uniprot

P15408

Synonyms:

Fos-related antigen 2 (FRA-2), FOSL2, FRA2


Immunogen:

Recombinant Human Fos-related antigen 2 protein (187-326AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43965 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung cancer using PACO43965 at dilution of 1:100.