

DHX58 Antibody



PACO43968

Product Information

Size:

50ul

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

DHX58

Uniprot

Q96C10

Synonyms:

Probable ATP-dependent RNA helicase DHX58 (EC 3.6.4.13) (Probable ATP-dependent helicase LGP2) (Protein D11Lgp2 homolog) (RIG-I-like receptor 3) (RLR-3) (RIG-I-like receptor LGP2) (RLR), DHX58, D11LGP2E LGP2


Immunogen:

Recombinant Human Probable ATP-dependent RNA helicase DHX58 protein (389-678AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: DHX58 antibody at 12 μ g/ml. Lane 1: 293T whole cell lysate. Lane 2: Rat kidney tissue. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 77 kDa. Observed band size: 77 kDa.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43968 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43968 at dilution of 1:100.