

CST2 Antibody



PACO44027

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Thiol protease inhibitor.

Gene ID:

CST2

Uniprot

P09228

Synonyms:

Cystatin-SA (Cystatin-2) (Cystatin-S5), CST2


Immunogen:

Recombinant Human Cystatin-SA protein (20-141AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded human salivary gland tissue using PACO44027 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44027 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO44027 at dilution of 1:100.