

TFF2 Antibody



PACO44042

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Inhibits gastrointestinal motility and gastric acid, secretion. Could function as a structural component of gastric mucus, possibly by stabilizing glycoproteins in the mucus gel through interactions with carbohydrate side chains.

Gene ID:

TFF2

Uniprot

Q03403

Synonyms:

Trefoil factor 2 (Spasmodysin) (Spasmolytic polypeptide) (SP), TFF2, SML1


Immunogen:

Recombinant Human Trefoil factor 2 protein (24-129AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO44042 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO44042 at dilution of 1:100.