

CNTNAP2 Antibody



PACO44046

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in the formation of functional distinct domains critical for saltatory conduction of nerve impulses in myelinated nerve fibers. Seems to demarcate the juxtaparanodal region of the axo-glia junction.

Gene ID:

CNTNAP2

Uniprot

Q9UHC6

Synonyms:

Contactin-associated protein-like 2 (Cell recognition molecule Caspr2), CNTNAP2, CASPR2 KIAA0868

Immunogen:

Recombinant Human Contactin-associated protein-like 2 protein (180-400AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A