

CRYBB1 Antibody



PACO44177

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Crystallins are the dominant structural components of the vertebrate eye lens.

Gene ID:

CRYBB1

Uniprot

P53674

Synonyms:

Beta-crystallin B1 (Beta-B1 crystallin), CRYBB1


Immunogen:

Recombinant Human Beta-crystallin B1 protein (1-252AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: CRYBB1 antibody at 1.58 μ g/ml. Lane 1: MCF-7 whole cell lysate. Lane 2: Mouse kidney tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 28 kDa. Observed band size: 28 kDa.