
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

extracellular exosome, extracellular region, extracellular space, fibrillar collagen trimer, Golgi lumen, lysosomal lumen, proteinaceous extracellular matrix, collagen binding, extracellular matrix structural constituent, collagen fibril organization.

Gene ID:

LUM

Uniprot

P51884

Synonyms:

Lumican (Keratan sulfate proteoglycan lumican) (KSPG lumican), LUM, LDC SLRR2D


Immunogen:

Recombinant Human Lumican protein (19-338AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO44182 at dilution of 1:100.