

TAF12 Antibody



PACO44204

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

TAFs are components of the transcription factor IID (TFIID) complex, PCAF histone acetylase complex and TBP-free TAFII complex (TFTC). TAFs components-TIIFD are essential for mediating regulation of RNA polymerase transcription.

Gene ID:

TAF12

Uniprot

Q16514

Synonyms:

Transcription initiation factor TFIID subunit 12 (Transcription initiation factor TFIID 20/15 kDa subunits) (TAFII-20/TAFII-15) (TAFII20/TAFII15), TAF12, TAF15 TAF2J TAFII20


Immunogen:

Recombinant Human Transcription initiation factor TFIID subunit 12 protein (1-161AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44204 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO44204 at dilution of 1:100.