

CDC7 Antibody



PACO44213

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Seems to phosphorylate critical substrates that regulate the G1/S phase transition and/or DNA replication. Can phosphorylates MCM2 and MCM3.

Gene ID:

CDC7

Uniprot

O00311

Synonyms:

Cell division cycle 7-related protein kinase (CDC7-related kinase) (HsCdc7) (huCdc7) (EC 2.7.11.1), CDC7, CDC7L1


Immunogen:

Recombinant Human Cell division cycle 7-related protein kinase protein (275-574AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44213 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44213 at dilution of 1:100.