

CRYGD Antibody



PACO44214

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Crystallins are the dominant structural components of the vertebrate eye lens.

Gene ID:

CRYGD

Uniprot

P07320

Synonyms:

Gamma-crystallin D (Gamma-D-crystallin) (Gamma-crystallin 4), CRYGD, CRYG4

Immunogen:

Recombinant Human Gamma-crystallin D protein (1-174AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A