

SERPIND1 Antibody



PACO44220

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Thrombin inhibitor activated by the glycosaminoglycans, heparin or dermatan sulfate. In the presence of the latter, HC-II becomes the predominant thrombin inhibitor in place of antithrombin III (AT-III). Also inhibits chymotrypsin, but in a glycosaminoglycan-independent manner.

Gene ID:

SERPIND1

Uniprot

P05546

Synonyms:

Heparin cofactor 2 (Heparin cofactor II) (HC-II) (Protease inhibitor leuserpin-2) (HLS2) (Serp D1), SERPIND1, HCF2


Immunogen:

Recombinant Human Heparin cofactor 2 protein (20-290AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO44220 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human Plasma tissue using PACO44220 at dilution of 1:100.