

BNIP2 Antibody



PACO44230

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Implicated in the suppression of cell death. Interacts with the BCL-2 and adenovirus E1B 19 kDa proteins.

Gene ID:

BNIP2

Uniprot

Q12982

Synonyms:

BCL2/adenovirus E1B 19 kDa protein-interacting protein 2, BNIP2, NIP2


Immunogen:

Recombinant Human BCL2/adenovirus E1B 19 kDa protein-interacting protein 2 protein (125-280AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO44230 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO44230 at dilution of 1:100.