

PACO44232

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

TCIM

Uniprot

Q9NR00

Synonyms:

Transcriptional and immune response regulator (Thyroid cancer protein 1) (TC-1), TCIM, C8orf4 TC1


Immunogen:

Recombinant Human Transcriptional and immune response regulator protein (1-106AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO44232 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung cancer using PACO44232 at dilution of 1:100.