

PACO44241

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Isoform 1: May be involved in RNA splicing.

Gene ID:

NOL3

Uniprot

O60936

Synonyms:

Nucleolar protein 3 (Apoptosis repressor with CARD) (Muscle-enriched cytoplasmic protein) (Myp) (Nucleolar protein of 30 kDa) (Nop30), NOL3, ARC NOP


Immunogen:

Recombinant Human Nucleolar protein 3 protein (1-208AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44241 at dilution of 1:100.