
Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

Part of the MIS12 complex, which is required for kinetochore formation and spindle checkpoint activity. Required to target ZW10 to the kinetochore at prometaphase.

Gene ID:

ZWINT

Uniprot

O95229

Synonyms:

ZW10 interactor (ZW10-interacting protein 1) (Zwint-1), ZWINT


Immunogen:

Recombinant Human ZW10 interactor protein (1-277AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: ZWINT antibody at 8.93µg/ml. Lane 1: Hela whole cell lysate. Lane 2: Mouse liver tissue. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 32, 26 kDa. Observed band size: 32 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44251 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO44251 at dilution of 1:100.