

# ADAD1 Antibody



PACO44254

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Plays a role in spermatogenesis. Binds to RNA but not to DNA.

**Gene ID:**

ADAD1

**Uniprot**

Q96M93

**Synonyms:**

Adenosine deaminase domain-containing protein 1 (Testis nuclear RNA-binding protein), ADAD1, TENR

**Immunogen:**


Recombinant Human Adenosine deaminase domain-containing protein 1 protein (1-280AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human breast cancer using PACO44254 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44254 at dilution of 1:100.