

ANKRD36B Antibody



PACO44278

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ANKRD36B

Uniprot

Q8N2N9

Synonyms:

Ankyrin repeat domain-containing protein 36B (CLL-associated antigen KW-1), ANKRD36B

Immunogen:

Recombinant Human Ankyrin repeat domain-containing protein 36B protein (1-230AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A