

ELF1 Antibody



PACO44283

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcription factor that activates the LYN and BLK promoters. Appears to be required for the T-cell-receptor-mediated trans activation of HIV-2 gene expression. Binds specifically to two purine-rich motifs in the HIV-2 enhancer.

Gene ID:

ELF1

Uniprot

P32519

Synonyms:

ETS-related transcription factor Elf-1 (E74-like factor 1), ELF1


Immunogen:

Recombinant Human ETS-related transcription factor Elf-1 protein (1-210AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO44283 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44283 at dilution of 1:100.