

SKAP2 Antibody



PACO44304

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in B-cell and macrophage adhesion processes. In B-cells, may act by coupling the B-cell receptor (BCR) to integrin activation. May play a role in src signaling pathway.

Gene ID:

SKAP2

Uniprot

O75563

Synonyms:

Src kinase-associated phosphoprotein 2 (Pyk2/RAFTK-associated protein) (Retinoic acid, induced protein 70) (SKAP55 homolog) (SKAP-55HOM) (SKAP-HOM) (Src family-associated phosphoprotein 2) (Src kinase-associated phosphoprotein 55-related protein) (Src-associated adapter protein with PH and SH3 domains), SKAP2, PRAP RA70 SAPS SCAP2 SKAP55R


Immunogen:

Recombinant Human Src kinase-associated phosphoprotein 2 protein (1-120AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44304 at dilution of 1:100.