
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Plays a role in the control of cell shape and motility in the trabecular meshwork.

Gene ID:

COCH

Uniprot

O43405

Synonyms:

Cochlin (COCH-5B2), COCH, COCH5B2


Immunogen:

Recombinant Human Cochlin protein (20-260AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO44325 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO44325 at dilution of 1:100.