

DLX2 Antibody



PACO44331

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis.

Gene ID:

DLX2

Uniprot

Q07687

Synonyms:

Homeobox protein DLX-2, DLX2


Immunogen:

Recombinant Human Homeobox protein DLX-2 protein (1-140AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO44331 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO44331 at dilution of 1:100.