

PPP2R1B Antibody



PACO44408

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:100-1:500

Protein Background:

The PR65 subunit of protein phosphatase 2A serves as a scaffolding molecule to coordinate the assembly of the catalytic subunit and a variable regulatory B subunit.

Gene ID:

PPP2R1B

Uniprot

P30154

Synonyms:

Serine/threonine-protein phosphatase 2A 65 kDa regulatory subunit A beta isoform (PP2A subunit A isoform PR65-beta) (PP2A subunit A isoform R1-beta), PPP2R1B

Immunogen:

Recombinant Human Serine/threonine-protein phosphatase 2A 65 kDa regulatory subunit A β isoform protein (402-601AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: PPP2R1B antibody at 3.33 μ g/ml + Jurkat whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 67, 74, 62, 53 kDa. Observed band size: 67 kDa.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO44408 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44408 at dilution of 1:100.