

VPS35 Antibody



PACO44422

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

VPS35

Uniprot

Q96QK1

Synonyms:

Vacuolar protein sorting-associated protein 35 (hVPS35) (Maternal-embryonic 3) (Vesicle protein sorting 35), VPS35, MEM3


Immunogen:

Recombinant Human Vacuolar protein sorting-associated protein 35 protein (1-250AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO44422 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44422 at dilution of 1:100.