

RPS27 Antibody



PACO44454

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytosol, cytosolic small ribosomal subunit, nucleoplasm, nucleus, ribosome, DNA binding, poly(A) RNA binding, structural constituent of ribosome, zinc ion binding, cell proliferation.

Gene ID:

RPS27

Uniprot

P42677

Synonyms:

40S ribosomal protein S27 (Metallopan-stimulin 1) (MPS-1) (Small ribosomal subunit protein eS27), RPS27, MPS1


Immunogen:

Recombinant Human 40S ribosomal protein S27 protein (1-84AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO44454 at dilution of 1:100.