
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in regulation of vitamin K-dependent carboxylation of multiple N-terminal glutamate residues. Seems to inhibit gamma-carboxylase GGCX. Binds 7 calcium ions with a low affinity.

Gene ID:

CALU

Uniprot

O43852

Synonyms:

Calumenin (Crocabin) (IEF SSP 9302), CALU


Immunogen:

Recombinant Human Calumenin protein (20-315AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO44461 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO44461 at dilution of 1:100.