

CD200 Antibody



PACO44462

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Costimulates T-cell proliferation. May regulate myeloid cell activity in a variety of tissues.

Gene ID:

CD200

Uniprot

P41217

Synonyms:

OX-2 membrane glycoprotein (CD antigen CD200), CD200, MOX1 MOX2


Immunogen:

Recombinant Human OX-2 membrane glycoprotein protein (30-240AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO44462 at dilution of 1:100.