

CD200 Antibody



PACO44463

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Costimulates T-cell proliferation. May regulate myeloid cell activity in a variety of tissues.

Gene ID:

CD200

Uniprot

P41217

Synonyms:

OX-2 membrane glycoprotein (CD antigen CD200), CD200, MOX1 MOX2

Immunogen:

Recombinant Human OX-2 membrane glycoprotein protein (30-240AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: CD200antibody at 0.71 μ g/ml + Mouse kidney tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 32, 31, 33 kDa. Observed band size: 32 kDa.