

TMED10 Antibody



PACO44471

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

TMED10

Uniprot

P49755

Synonyms:

Transmembrane emp24 domain-containing protein 10 (21 kDa transmembrane-trafficking protein) (S31III125) (S31I125) (Tmp-21-I) (Transmembrane protein Tmp21) (p23) (p24 family protein delta-1) (p24delta1) (p24delta), TMED10, TMP21


Immunogen:

Recombinant Human Transmembrane emp24 domain-containing protein 10 protein (35-190AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: TMED10 antibody at 4.88 μ g/ml. Lane 1: HeLa whole cell lysate. Lane 2: Mouse liver tissue. Lane 3: Mouse lung tissue. Lane 4: 293T whole cell lysate. Lane 5: NIH/3T3 whole cell lysate. Lane 6: U87 whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 25 kDa. Observed band size: 25 kDa.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44471 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO44471 at dilution of 1:100.