

PACO44476

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Supplies the nucleotide substrate for thymidylate synthetase.

**Gene ID:**

DCTD

**Uniprot**

P32321

**Synonyms:**

Deoxycytidylate deaminase (EC 3.5.4.12) (dCMP deaminase), DCTD

**Immunogen:**


Recombinant Human Deoxycytidylate deaminase protein (1-178AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver cancer using PACO44476 at dilution of 1:100.