

BTNL2 Antibody



PACO44487

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

Negative regulator of T-cell proliferation.

Gene ID:

BTNL2

Uniprot

Q9UIR0

Synonyms:

Butyrophilin-like protein 2 (BTL-II), BTNL2


Immunogen:

Recombinant Human Butyrophilin-like protein 2 protein (246-455AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: BTNL2 antibody at 4.12 μ g/ml. Lane 1: Mouse small intestine tissue. Lane 2: K562 whole cell lysate. Lane 3: HepG2 whole cell lysate. Lane 4: HL60 whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 51, 44, 31, 28, 20, 29 kDa. Observed band size: 51 kDa.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO44487 at dilution of 1:100.