

EIF2B2 Antibody



PACO44495

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP.

Gene ID:

EIF2B2

Uniprot

P49770

Synonyms:

Translation initiation factor eIF-2B subunit beta (S20I15) (S20III15) (eIF-2B GDP-GTP exchange factor subunit beta), EIF2B2, EIF2BB


Immunogen:

Recombinant Human Translation initiation factor eIF-2B subunit β protein (1-351AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44495 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO44495 at dilution of 1:100.