

PSMD8 Antibody



PACO44520

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Acts as a regulatory subunit of the 26S proteasome which is involved in the ATP-dependent degradation of ubiquitinated proteins. Necessary for activation of the CDC28 kinase.

Gene ID:

PSMD8

Uniprot

P48556

Synonyms:

26S proteasome non-ATPase regulatory subunit 8 (26S proteasome regulatory subunit RPN12) (26S proteasome regulatory subunit S14) (p31), PSMD8


Immunogen:

Recombinant Human 26S proteasome non-ATPase regulatory subunit 8 protein (101-350AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO44520 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO44520 at dilution of 1:100.